



Employing Barcode High-Resolution Melting Technique for Authentication of Apricot Cultivars

Kaan HÜRKAN^{a,b}

^aDepartment of Agricultural Biotechnology, Faculty of Agriculture, Iğdır University, 76000, TURKEY

^bIğdır University, Research Laboratory Practice and Research Center, 76000, Iğdır, TURKEY

ARTICLE INFO

Research Article

Corresponding Author: Kaan HÜRKAN, E-mail: kaan.hurkan@igdir.edu.tr, kaanhurkan@hotmail.com

Received: 29 January 2021 / Revised: 11 May 2021 / Accepted: 12 May 2021 / Online: 25 March 2022

ABSTRACT

Fast, accurate and affordable identification of food products is important to ensure authenticity and safety. There are various apricot (*Prunus armeniaca* L.) cultivars are being produced in Turkey. Each cultivar differs in quality and purpose of use. In this paper, we aimed to develop an easy and reliable method, Barcode High-Resolution Melting (Bar-HRM), to distinguish apricot cultivars. We designed and tested novel Bar-

HRM primer sets HRM-ITS1 and HRM-ITS2, targeting the most popular barcoding region ITS1 and ITS2, specific to apricot cultivars. According to the results, HRM analysis distinguished 31 cultivars of 35 for ITS1, and 35 for ITS2. We recommend using ITS2 barcode region, amplified with using HRM-ITS2 primer set, for Bar-HRM analysis of different apricot cultivars.

Keywords: *Prunus armeniaca*, Internal Transcribed Spacer, Bar-HRM, Molecular Identification

© Ankara University, Faculty of Agriculture

1. Introduction

Apricot (*Prunus armeniaca* L.) is an important drupe fruit with a rich gene pool that makes the plant capable to adapt from Siberia to North Africa and Greece to China (Mehlenbacher et al. 1991). Turkey leads to world apricot production with 677,000 tons on average (FAOSTAT 2020). Apricot production is specialised especially in Malatya, Erzincan and Iğdır regions in Turkey (Ercisli 2004). Because apricot has been cultivated in Inner Anatolia since ancient times, some cultivars adapted to different climatic areas of the country i.e. Hacıhaliloğlu in Malatya, Hasanbey in Erzincan and Şalak in Iğdır (Güleryüz et al. 1997). According to the Turkish Apricot Research Institute data, there are 28 registered cultivars and numerous genotypes in Turkey including the Protected Designation of Origin (PDO) cultivars Hacıhaliloğlu, Hasanbey, Kabaş, Soğancı, Çataloğlu, Çöloğlu, Şalak (Anonymous 2021).

The economic importance of the PDO products is higher than the common products. Both producers and consumers must ensure PDO authenticity against adulteration. This makes “tracing the original food product” is crucial.

DNA barcoding, a method based on comparing nucleotide sequences of a specific DNA fragment, is a widely used tool for identifying the species, reconstructing the phylogenetics, assessing the biodiversity since the early 2000s (Cheng et al. 2016). Different DNA regions have been using as DNA barcodes for various plant groups such as *rbcL*, *matK*, *ycf1*, *LEAFY* and the most used Internal Transcribed Spacer (ITS) region (Won & Renner 2005; Chase et al. 2007; Kress & Erickson 2007; Hollingsworth 2011; Li et al. 2011; Dong et al. 2015). It is also confirmed that ITS barcode region distinguishes apricot cultivars successfully (Hürkan 2020). Although the DNA barcoding technique is very useful for species identification, it prolongs the workflow, increases the expenses due to post-PCR sequencing analysis, and needs a bioinformatics background for the researcher. Recently, DNA barcoding technique supported with other techniques to overcome these disadvantages. In this paper, we urged on the High Resolution Melting (HRM) technique coupled with DNA barcoding region to develop a cost-effective method to discriminate closely related apricot cultivars.

High Resolution Melting analysis is a technique that used for genotyping by discriminating DNA sequence differences of Single Nucleotide Polymorphisms (SNPs) and sequence length polymorphisms in PCR products (Zhou et al. 2005). In this technique, the shape of the melting transitions of the PCR products is being acquiesced continuously. HRM is a more powerful, cheaper and easier technique than other approaches since requires neither post-PCR processing nor bioinformatics skills. Moreover, it is faster and more economical since being a sequencing-free method. Although HRM technique has majorities over conventional methods, it had some disadvantages that could directly influence the results. During the early stages of the HRM,

SYBR® Green was used as dsDNA binding dye, which inhibits the DNA polymerase. Therefore, SYBR™ Green dye does not allow discriminating closely related genotypes that have small sequence variations on HRM (Reed et al. 2007). However, saturation dyes e.g. LC Green™, SYTO9™ or Eva Green™, which is selected for this study, do not affect the DNA polymerase performance even at high concentrations (Vossen et al. 2009). The discriminating power of HRM directly related to the markers used in the analysis. In recent researches, there are two marker types are being used for HRM analysis, microsatellites and DNA barcoding regions. Both have advantages and disadvantages. While microsatellites are usually organism-group specified and need more experience to design and mine the primers, DNA barcoding regions are universal.

Combining the discrimination success of DNA barcoding and the easiness of HRM results in successful identification of food products such as PDO cheese (Ganopoulos et al. 2013), codfish species (Shi et al. 2020), poisonous plants (Thongkhao et al. 2020), medicinal plants (Li et al. 2016; Sun et al. 2017; Mishra et al. 2018), and bean crops (Madesis et al. 2012). In a recent study, 16 SSR markers from the literature were used to characterise the wild apricot genotypes from Nevşehir region (Turkey) (Bakır et al. 2019).

Because of the advantages above, herein, we developed and tested a Bar-HRM based method, using specific primer sets targeting the most used barcoding region ITS1 and ITS2, to rapid, affordable, reliable and quantitative identification of apricot cultivars and apricot products.

2. Material and Methods

2.1. Samples collection and DNA extraction

The Republic of Turkey, Ministry of Agriculture and Forestry Apricot Research Institute (Malatya, Turkey) kindly provided fresh leaf samples of the 35 apricot cultivars available in April 2020 (Table 1).

Table 1- Samples information, Cycle Threshold (C_q) values, average Melting Temperature (T_m) and genotyping confidence percentages

Sample Code	Colour	Cultivar Name	C _q ±SD		Average T _m (°C)±SD		Confidence (%) ²	
			ITS1	ITS2	ITS1	ITS2	ITS1	ITS2
AS	—	Aprikoz (Şalak)	9.49±0.21	11.66±0.12	87.70±0.06	90.56±0.04	100.00 ³	100.00 ³
AY	—	Alyanak ¹	9.55±0.18	11.47±0.85	87.34±0.06	90.20±0.15	69.06	63.17
AZ	—	Aprikoz ¹	12.77±0.04	14.70±0.55	87.60±0.06	90.26±0.08	93.91	85.67
CG	—	Çekirge 52	9.21±1.02	11.72±0.45	87.40±0.10	90.26±0.06	91.66	81.27
CT	—	Çataloğlu ¹	9.08±0.71	10.95±0.38	87.56±0.10	90.10±0.10	86.42	63.11
EA	—	Erken ağerik	9.99±0.11	11.53±0.65	87.36±0.06	90.14±0.06	66.64	68.92
EB	—	Ethembey ¹	10.00±0.07	12.12±0.47	87.34±0.10	90.16±0.12	74.15	62.24
HB	—	Hasanbey ¹	9.76±0.065	11.22±0.55	87.40±0.10	90.04±0.04	96.69 ³	59.31
HG	—	Hungarish Beste	8.57±1.07	10.70±0.48	87.50±0.06	90.26±0.06	91.67	78.07
HH	—	Hacıhaliloğlu ¹	8.98±0.66	10.62±0.11	87.50±0.06	90.16±0.10	85.92	76.57
HK	—	Hacıkız	10.27±0.41	12.22±0.09	87.50±0.10	90.10±0.18	94.19	63.78
HR	—	Hırmanlı Dırağı	9.82±0.94	12.18±0.16	87.40±0.06	90.10±0.05	96.54 ³	75.19
IA	—	İsmailağa	9.49±0.37	11.23±1.14	87.66±0.06	90.04±0.06	85.07	52.55
IB	—	İri bitirgen ¹	9.40±0.98	11.27±0.87	87.44±0.10	90.14±0.06	86.14	56.94
IH	—	İmrahor ¹	9.60±1.22	11.68±0.61	87.26±0.10	90.20±0.10	21.91	40.45
KA	—	Kabaaşı ¹	9.52±1.18	11.97±1.55	87.50±0.10	90.10±0.08	98.64 ³	80.55
KB	—	Karacabey ¹	9.15±0.57	11.36±0.78	87.60±0.06	90.26±0.10	71.48	73.32
KK	—	Kuru Kabuk ¹	9.93±0.098	10.93±0.08	87.50±0.10	90.24±0.11	20.54	30.47
PB	—	Precose de Boulbon	9.76±0.19	11.65±0.41	87.36±0.11	90.04±0.08	27.56	11.76
PL	—	Polenais	10.81±0.48	13.16±0.68	87.96±0.08	90.16±0.06	81.46	78.88
PM	—	Paşa Mişmişi	10.92±0.98	12.58±0.49	87.44±0.14	90.10±0.10	58.15	81.00
PV	—	Pavlot	9.58±1.54	11.71±1.24	87.56±0.08	90.30±0.14	78.72	88.71
RK	—	Rakowsky	10.46±0.64	12.98±0.34	87.46±0.06	90.24±0.10	48.14	45.63
RY	—	Royal	9.82±0.15	11.86±1.23	87.60±0.14	90.24±0.09	77.44	73.40
SC	—	Soğancı ¹	9.99±0.45	12.10±0.24	87.90±0.10	90.20±0.18	63.38	82.10
SE	—	Stark Early Orange ¹	10.01±0.11	12.40±0.42	88.06±0.06	90.36±0.04	78.00	81.35
SK	—	Şekerpare ¹	10.48±1.02	12.95±0.35	87.64±0.12	90.50±0.10	77.76	75.21
SM	—	Şam ¹	10.52±0.12	12.67±0.57	87.70±0.02	90.46±0.06	86.07	94.84
TE	—	Tokalolu (Konya Ereğli)	10.12±0.98	12.54±1.74	87.80±0.13	90.44±0.06	49.92	70.13
TF	—	Turfanda (Eski Malatya) ¹	10.32±0.41	12.52±0.35	87.84±0.17	90.36±0.08	45.72	87.63
TK	—	Tekeler	10.16±0.48	12.70±0.54	87.90±0.09	90.34±0.10	29.86	82.66
TO	—	Tokalolu (Erzincan)	12.07±0.94	14.58±0.38	87.90±0.06	90.54±0.12	45.60	39.47
TY	—	Tokalolu (Yalova)	10.75±0.89	13.06±0.68	88.00±0.14	90.84±0.06	44.10	59.89
WD	—	Wilson Delicouse	10.93±0.87	13.90±1.56	88.16±0.05	91.10±0.06	1.35	5.90
YN	—	Yeğen (Eski Malatya)	11.23±0.14	14.18±0.34	88.40±0.06	91.30±0.14	0.30	1.06

¹: Cultivars are registered to the Republic of Turkey, Ministry of Agriculture and Forestry Apricot Research Institute (Malatya, Turkey); ²: Confidence percentage calculated regarding AS; ³: Confidence value >95% infers no genotyping success, <95% successful genotyping

To extract DNA, we used ~100 mg of fresh leaf tissues and followed the modified CTAB protocol described in the literature (Aydın et al. 2018). The DNA concentration was measured by NanoDrop (Maestrogen) and integrity was confirmed on agarose gel electrophoresis. We normalised the concentration of all the DNA samples to 10 ng μL^{-1} . The DNA samples were stored at -20 °C for further analyses.

2.2. Primers mining

Since there is no available nucleotide sequence of the studied cultivars available on GenBank, we designed the primer sets for ITS1 and ITS2 barcode regions according to the GenBank apricot nucleotide sequence records with the accession number MT072696, EF211085, EF211084, EF211083 and MG735482. We downloaded the GenBank formatted files with annotations, imported them to the Geneious R8 software (Kearse et al. 2012) and evaluated them for quality and variable characters. Sequence characteristics were also analysed using MEGA Version X (Kumar et al. 2018). Then, the sequences were aligned (Geneious Alignment Tool) and the novel primer sets were designed by considering the variable positions, GC content, expected amplicon size and melting temperature (T_m) in the same software (Table 2). We in silico confirmed the specificity of the primers on the Primer-BLAST tool in the National Center for Biotechnology Information (NCBI).

Table 2- Primers information

Primer Name	Target Region	Sequence (5'→3')	T_m (°C)	GC Content (%)	Expected (Detected) Amplicon Size (bp)
HRM-ITS1	ITS1	F: TCGAAACCTGCCTAGCAGAA	59.0	50.0	121 (110)
		R: CGCGCTCTCTCGTTCAAGTT	61.0	55.0	
HRM-ITS2	ITS2	F: GAAGCCATTAGGCCGAGGG	60.2	63.2	144 (136)
		R: GGTTTCGCAACCACCGATTG	60.4	55.0	

Detected amplicon sizes by sequencing were given in parenthesis; T_m : Melting temperature; bp: Base pair

2.3. HRM analysis and sequencing of PCR products

We performed the HRM amplifications using Rotor-Gene-Q 5plex thermal cycler (Qiagen, USA) with a 72-well carousel. The HRM reaction mix was prepared as 5 μL Luminaris Colour HRM Master Mix (Thermo Scientific, USA), 0.5 μL of 10 mM each primer (Sentebiolab, Turkey), 10 ng template and nuclease-free water to 10 μL total volume. The cycling protocol was 95 °C 10 min initial denaturation followed by 40 cycles of 95 °C 10 s denaturation, 60 °C 30 s annealing, and 72 °C 30 s extension. Data acquired following each extension step. We added 95 °C 30 s and 50 °C 30 s steps for heteroduplex formation to the end of the cycle. We performed HRM immediately after the amplification in increments of 0.1 °C s^{-1} from 75 °C to 95 °C and data acquired continuously. All the reactions were performed as triplicates and no template control (NTC) was included in the reactions.

For HRM data analysis, we used the Rotor-Gene Q 2.3.5 (Qiagen, Germany) software. We calculated the Cycle Threshold (C_q) values by comparative quantification method of the software, Melting Temperatures (T_m), and normalised the HRM curves by removing the background fluorescence. The difference plots were generated regarding AS. Then, the software calculated the Genotype Confidence Percentage (GCP) for each cultivar. We set the confidence threshold to 95%.

To validate HRM analysis results, four ITS1 PCR products for the cultivars AS, HB, HR, KA which could not be separated by HRM (confidence percentage >95% threshold) and five ITS2 PCR products for the cultivars AS, AZ, PV, SM, TF which were separated by HRM (confidence percentage >85%). We sent the PCR products to the direct sequencing (Macrogen Europe) were performed in both directions using the same primers used for HRM amplifications by on ABI 3730xl System. We used Geneious R8 (Kearse et al. 2012) software for assembling the sequences and generating the consensus sequences. The sequences were aligned by the Geneious Aligner algorithm in the same software.

3. Results

3.1. DNA extraction and data mining

We obtained a sufficient amount of DNA from the samples with the followed extraction protocol. The DNA concentrations ranged from 29.69 to 187.29 ng μL^{-1} , and the A260/230 ratio ranged from 1.540 to 2.000.

We first analysed the sequences of ITS1 and ITS2 regions from the cultivars retrieved from GenBank to generate an initial comparison of the characteristics (Table 3). The length of the ITS1 region was identical among the cultivars as 215 bp. Ten variable sites (4.65%), 205 conserved sites (95.35%), and one parsimony informative sites were observed for the ITS1 region. ITS2 region was relatively longer (277.8±0.4 bp) and three folds by the variable sites (37 sites equals 13.32%) compared to ITS1. The conserved sites were lower 241 (86.75%) than ITS1, and there was no parsimony-informative site on ITS2. The average GC contents were almost the same for both regions (Figure 1). We would like to include the coding 5.8S region in the comparison

table to demonstrate why this coding region is not suitable for any genotyping analysis for apricot since all the sequence was conserved and there was no variable site available.




Figure 1- Comparison of average length (A), variable sites (B) and average GC per cent (C) of apricot cultivar ITS1, 5.8S and ITS2 sequences retrieved from GenBank.

Table 3- Characteristics of sequences retrieved from GenBank and sequencing of the HRM-PCR products

<i>Statistics</i>	<i>GenBank</i>			<i>HRM Amplicons</i>	
	<i>ITS1</i>	<i>5.8S</i>	<i>ITS2</i>	<i>ITS1</i>	<i>ITS2</i>
Average ungapped length (bp)+SD	215±0	103±0	277.8±0.4	110±0	136±0
Variable sites	10	0	37	0	6
Conserved sites	205	103	241	110	130
Parsimony informative sites	1	0	0	0	1
Average GC content (%)	63.4	47.6	64.4	60.9	61.8

3.2. HRM and sequencing results

ITS1 primer set amplified the expected size products, approximately 160 bp long. The software calculated the threshold cycles (C_q) and melting temperatures (T_m) of the amplicons (Table 1). The C_q values for ITS1 ranged from 8.57 to 12.77, and this marker generated 31 different, 4 identical (AS-IA, HB-PB, HR-RY and EA-SC) C_q values. The T_m values of the amplicons ranged from 87.26 to 88.40, and the software calculated 20 different T_m values. Although only 20 different T_m calculated of the 35 cultivars, HRM analysis successfully distinguished 32 cultivars (Figure 2). The cultivars HB, HR and KA could not be distinguished, and the GCPs were above 95% (96.69, 96.54 and 98.64, respectively). Figure 1 displays the HRM results utilizing difference and normalised plots. All the cultivars, except HB, HR and KA, yielded unique HRM profiles. So that, the software could distinguish the 32 cultivars by GCP calculations.




Figure 2- HRM analysis result of the apricot cultivars with the ITS1 nuclear marker. A) Normalised difference plot of 35 cultivars. B) Assigned genotypes plot drawn as AS control. Cut off value for Genotype confidence is 95%. Colour codes of the genotypes were indicated in Table 1. Normalisation ranges were adjusted to 84 °C (leading) and 90 °C (trailing).

ITS2 marker yielded better results than the ITS1. ITS2 primer set amplified approximately 170 bp long PCR products, as expected. The C_q values were ranged from 10.62 to 14.70, and each value was unique (Table 1). The T_m values ranged from 90.04 to 91.30, and the amplicons had 18 unique T_m values. HRM analysis of ITS2 distinguished all the cultivars by melting shapes (Figure 3), as we aimed. In contrast to conventional T_m difference analysis, the normalisation of the fluorescence by the software distinguished all the cultivars by melting curve shapes for ITS2.




Figure 3- HRM analysis result of the apricot cultivars with the ITS2 nuclear marker. A) Normalised difference plot of 35 cultivars. B) Assigned genotypes plot drawn as AS control. Cut off value for Genotype confidence is 95%. Colour codes of the genotypes were indicated in Table 1. Normalisation ranges were adjusted to 86.6 °C (leading) and 92.2 °C (trailing).

We analysed the amplicon sequences to validate the results of the HRM analysis. After trimming the primer binding sites, we obtained 110 bp sequence for ITS1 and 136 bp sequence for the ITS2 region (Table 2). Sequencing results showed while ITS1 sequences of the AS, HB, HR, KA were identical, six variable sites were detected for ITS2 sequences of AS, AZ, PV, SM, TF (Figure 4, Figure 5 and Table 3).




Figure 4- Alignment of the ITS1 target region for Bar-HRM analysis. The amplicon is 110 bp long. There are no nucleotide variable sites. The numbers indicate the position of the amplicon according to the complete ITS region.




Figure 5- Alignment of the ITS2 target region for Bar-HRM analysis. The amplicon is 136 bp long. Nucleotide variable sites are indicated with colours. The numbers indicate the position of the amplicon according to the complete ITS region.

4. Discussion

The ITS region is a suitable DNA barcoding region for plants (Kress et al. 2005), but this region has some problems such as duplications, paralogous copies and causes pseudogenes on some plant groups (Chase et al. 2007). The success of the ITS as a barcoding region and an HRM marker was reported for apricot cultivars (Hürkan 2020), Fabaceae (Gao & Chen 2009), *Artemisia* spp. (Song et al. 2016), and *Medicago lupulina* and *Trifolium pratense* (Ganopoulos et al. 2012). The ITS region consists of two non-coding, ITS1 and ITS2, and one coding, 5.8S, parts (Cheng et al. 2016). Systematic researches have proposed the ITS2 could be the core DNA barcode due to the region has high interspecific divergence (Xin et al. 2013). In this study, the basic comparison of the ITS parts showed parallel results to the literature as ITS2 was the most variable (13.22%) region followed by ITS1 (4.65%) for the studied apricot cultivars. Although the 5.8S gene coding region was well characterised for inter-specific level, the region has very limited nucleotide variations in deeper levels since it is a conserved region (Hershkovitz & Lewis 1996). Supporting the literature, the 5.8S region showed no variable regions for the studied apricot cultivars.

The optimal amplicon length for HRM analysis should be shorter than 300 bp (Reed & Wittwer 2004). Shorter amplicons emphasize Single Nucleotide Polymorphisms (SNPs) in HRM analysis. Both ITS1 and ITS2 primers worked fine for each studied cultivar and yielded approximately 160 and 170 bp amplicons, respectively, which were in the “ideal” range for HRM analysis. The melting temperature of the PCR products depends on both sequence length and nucleotide content. However, solely T_m values are not reliable for discrimination the organisms as seen in the results. ITS1 primer yielded 20 unique T_m values, and ITS2 was only 18. During the HRM reaction, following the PCR amplification, the thermal cycler applies temperature increment to PCR products to denature the double-strand DNA, while a detector continuously tracking the fluorescence change. Thus, the software considers not only T_m values but also the melting shapes of the amplicons. This provides the software with better discrimination ability (Reja et al. 2010, Reed & Wittwer 2004). After the normalisation of the melting shapes on the software, HRM analysis distinguished 32 cultivars of the 35 for ITS1, and 35 for ITS2. This result supports the comparison of the variable sites (Table 3 and Figure 1), and the literature (Ganopoulos et al. 2012; Song et al. 2016; Pereira et al. 2018; Mishra et al. 2018). ITS2 region had three times more variable sites than the ITS1 and yielded better discrimination result.

The sequencing results were parallel to both comparisons of the sequences retrieved from GenBank (Table 3) and the HRM analysis results. Although obtained amplicon lengths were shorter than expected, we found the variable sites for ITS2, which was necessary to distinguish the cultivars. Sequencing results showed ITS2 (six variable sites) was the more variable region than ITS1 (no variable sites). The lack of variable sites on ITS1 sequences for the cultivars, which have a confidence percentage higher than 95%, clearly explains why HRM analysis could not distinguish the cultivars. The identical sequences resulted in similar melting profiles. In contrast to ITS1, ITS2 had six variable sites on the sequences. All these variations we detected were transversion mutations. There was no insertion or deletion. The selected ITS2 samples for sequencing were the samples that have more than 85% confidence percentage. The sequencing result of ITS2 validated the sensitivity and reliability of the Bar-HRM analysis by distinguishing genetically closely related cultivars.

Standard DNA barcoding based Sanger-sequencing is relatively expensive. In this study, we used a commercial HRM kit and a single HRM reaction was costed 0.31 USD, while single direction Sanger-sequencing was costed 5.5 USD. Next-generation sequencing (NGS) based genotyping is a comparable level on cost. However, it needs sophisticated workflow e.g. pre-sequencing library preparation, high-level computing for post-sequencing assembly of the reads. It also needs experienced and expert researcher on bioinformatics. Therefore, we believe a well-designed Bar-HRM assay would be a fast, robust and cost-effective way to genotype organism groups.

5. Conclusions

Bar-HRM is a cost-effective, fast and robust identification method. Moreover, no specialisation is needed for handling the data. In this study, ITS based Bar-HRM has proven for the identification of the 35 apricot cultivars. ITS2 primers set had the highest discrimination rate and can be used for the identification of various apricot cultivars. ITS2 marker can also be used for identification of *Prunus* species, authentication of apricot products.

Acknowledgements

We thank Cemil Ernim (Administrative and Technical Coordinator of Republic of Turkey Ministry of Agriculture and Forestry Apricot Research Institute) for sharing the apricot collection and anonymous referees for the improvement of the manuscript.

Funding

No funding available.

References

Anonymous (2021). Republic of Turkey Ministry of Agriculture and Forestry

- Apricot Research Institute (Tarım ve Orman Bakanlığı Kayısı Araştırma Enstitüsü Müdürlüğü). Retrieved in January, 15, 2021 from <https://arastirma.tarimorman.gov.tr/kayisi/Menu/6/Tescilli-Cesitlerimiz>
- Bakır M, Dumanoglu H, Erdoğlan V, Ernim C, Macit T (2019). Characterization of Wild Apricot (*Prunus armeniaca* L.) Genotypes Selected from Cappadocia Region (Nevşehir-Turkey) by SSR Markers. *Journal of Agricultural Sciences* 25(4): 498-507. DOI: 10.15832/ankutbd.457850
- Chase M W, Cowan R S, Hollingsworth P M & van den Berg C (2007). A proposal for a standardised protocol to barcode all land plants. *Taxon* 56(2): 295-299 <https://doi.org/10.1002/tax.562004>
- Cheng T, Xu C, Lei L, Li C, Zhang Y & Zhou S (2016). Barcoding the kingdom Plantae: new PCR primers for ITS regions of plants with improved universality and specificity. *Molecular Ecology Resources* 16(1): 138-149 <https://doi.org/10.1111/1755-0998.12438>
- Dong W, Xu C, Li C et al (2015). *ycf1*, the most promising plastid DNA barcode of land plants. *Scientific Reports* 5: 8348 <https://doi.org/10.1038/srep08348>
- Ercisli S (2004). A short review of the fruit germplasm resources of Turkey. *Genetic Resources and Accessed on Crop Evolution* 51: 419-435 <https://doi.org/10.1023/b:gres.0000023458.60138.79>
- FAOSTAT (2020). Apricot production. Retrieved in December, 12, 2020 from <http://www.fao.org/faostat/en/>
- Ganopoulos I, Madesis P & Tsaftaris A (2012). Universal ITS2 barcoding DNA region coupled with high-resolution melting (HRM) analysis for seed authentication and adulteration testing in leguminous forage and pasture species. *Plant Molecular Biology Reporter* 30(6): 1322-1328 <https://doi.org/10.1007/s11105-012-0453-3>
- Ganopoulos I, Sakaridis I, Argiriou A, Madesis P & Tsaftaris A (2013). A novel closed-tube method based on high resolution melting (HRM) analysis for authenticity testing and quantitative detection in Greek PDO Feta cheese. *Food chemistry* 141(2): 835-840 <https://doi.org/10.1016/j.foodchem.2013.02.130>
- Gao T & Chen S L (2009). Authentication of the medicinal plants in Fabaceae by DNA barcoding technique. *Planta Med* 75: 417-417 <https://doi.org/10.1055/s-2009-1216451>
- Gülyeryüz M, Ercisli S & Esitken A (1997). A study on characteristic features of apricot grown in Erzincan, Malatya and Iğdir provinces. In XI International Symposium on Apricot Culture 488: 165-170 <https://doi.org/10.17660/actahortic.1999.488.22>
- Hollingsworth P M, Graham S W & Little D P (2011). Choosing and using a plant DNA barcode. *PloS one* 6(5): 19254 <https://doi.org/10.1371/journal.pone.0019254>
- Hürkan K (2020). Analysis of Various DNA Barcodes on the Turkish Protected Designation of Origin Apricot "İğdir Kayısı" (*Prunus armeniaca* cv. Şalak). *Turkish Journal of Agriculture-Food Science and Technology* 8(9): 1982-1987 <https://doi.org/10.24925/turjaf.v8i9.1982-1987.3594>
- Kress W J & Erickson D L (2007). A two-locus global DNA barcode for land plants: the coding *rbcl* gene complements the non-coding *trnH-psbA* spacer region. *PLoS one* 2(6): 508 <https://doi.org/10.1371/journal.pone.0000508>
- Kress W J, Wurdack K J, Zimmer E A, Weigt L A, Janzen D H (2005). Use of DNA barcodes to identify flowering plants. *Proceedings of the National Academy of Sciences of the United States of America* 102: 8369-8374 <https://doi.org/10.1073/pnas.0503123102>
- Kumar S, Stecher G, Li M, Knyaz C & Tamura K (2018). MEGA X: Molecular Evolutionary Genetics Analysis across computing platforms. *Molecular Biology and Evolution* 35: 1547-1549 <https://doi.org/10.1093/molbev/msy096>
- Li D Z, Gao L M, Li H T et al (2011). Comparative analysis of a large dataset indicates that internal transcribed spacer (ITS) should be incorporated into the core barcode for seed plants. *Proceedings of the National Academy of Sciences of the United States of America* 108: 19641-19646 <https://doi.org/10.1073/pnas.1104551108>
- Li J, Song M, Xiong C, Zhao B & Sun W (2016). Application of barcode high-resolution melting for rapid authentication of the medicinal plant *Psammosilene tunicoides*. *Biotechnology & Biotechnological Equipment* 30(4): 790-796 <https://doi.org/10.1080/13102818.2016.1181988>
- Hershkovitz M A & Lewis L A (1996). Deep-level diagnostic value of the rDNA-ITS region. *Molecular Biology and Evolution* 13(9): 1276-1295 doi.org/10.1093/oxfordjournals.molbev.a025693
- Madesis P, Ganopoulos I, Anagnostis A & Tsaftaris A (2012). The application of Bar-HRM (Barcode DNA-High Resolution Melting) analysis for authenticity testing and quantitative detection of bean crops (Leguminosae) without prior DNA purification. *Food Control* 25(2): 576-582 <https://doi.org/10.1016/j.foodcont.2011.11.034>
- Mehlenbacher S A, Cociu V & Hough F L (1991). Apricots (*Prunus*). *Genetic Resources of Temperate Fruit and Nut Crops* 290: 65-110 <https://doi.org/10.17660/actahortic.1991.290.3>
- Mishra P, Shukla A K & Sundaresan V (2018). Candidate DNA barcode tags combined with high resolution melting (Bar-HRM) curve analysis for authentication of *Senna alexandrina* Mill. with Validation in Crude Drugs. *Frontiers in plant science* 9: 283 <https://doi.org/10.3389/fpls.2018.00283>
- Pereira L, Gomes S, Barrias S, Fernandes J R & Martins L P (2018). Applying high-resolution melting (HRM) technology to olive oil and wine authenticity. *Food Research International* 103: 170-181 <https://doi.org/10.1016/j.foodres.2017.10.026>
- Reed G H, Kent J O, & Wittwer C T (2007). High-resolution DNA melting analysis for simple and efficient molecular diagnostics. *Pharmacogenomics* 8(6): 597-608 <https://doi.org/10.2217/14622416.8.6.597> <https://doi.org/10.2217/14622416.8.6.597>
- Reed G H & Wittwer C T (2004). Sensitivity and specificity of single-nucleotide polymorphism scanning by high-resolution melting analysis. *Clinical Chemistry* 50: 1748-1754 <https://doi.org/10.1373/clinchem.2003.029751>
- Reja V, Kwok A, Stone G, Yang L, Missel A & Menzel C (2010). ScreenClust: Advanced statistical software for supervised and unsupervised high resolution melting (HRM) analysis. *Methods* 50: 10-14 <https://doi.org/10.1016/j.ymeth.2010.02.006>
- Shi R, Xiong X, Huang M, Xu W, Li Y, Cao M & Xiong X (2020). High resolution melting (HRM) analysis of a 12S rRNA mini barcode as a novel approach for codfish species authentication in processed fish products. *European Food Research and Technology* 1-9. <https://doi.org/10.1007/s00217-020-03456-5>
- Song M, Li J, Xiong C, Liu H & Liang J (2016). Applying high-resolution melting (HRM) technology to identify five commonly used *Artemisia* species. *Scientific reports* 6(1): 1-7 <https://doi.org/10.1038/srep34133>
- Sun W, Yan S, Li J, Xiong C, Shi Y, Wu L & Chen S (2017). Study of commercially available *Lobelia chinensis* products using Bar-HRM technology. *Frontiers in plant science* 8: 351 <https://doi.org/10.3389/fpls.2017.00351>
- Thongkhao K, Tungphatthong C, Phadungcharoen T & Sukrong S (2020). The use of plant DNA barcoding coupled with HRM analysis to differentiate edible vegetables from poisonous plants for food safety. *Food Control* 109: 106896 <https://doi.org/10.1016/j.foodcont.2019.106896>

- Vossen R H, Aten E, Roos A, & den Dunnen J T (2009). High-Resolution Melting Analysis (HRMA)—More than just sequence variant screening. *Human mutation*, 30(6): 860-866 <https://doi.org/10.1002/humu.21019>
- Won H & Renner S S (2005). The internal transcribed spacer of nuclear ribosomal DNA in the gymnosperm *Gnetum*. *Molecular Phylogenetics and Evolution* 36(3): 581-597 <https://doi.org/10.1016/j.ympev.2005.03.011>
- Xin T Y et al (2013). Super food *Lycium barbarum* (Solanaceae) traceability via an internal transcribed spacer 2 barcode. *Food Research International* 54: 1699-1704 <https://doi.org/10.1016/j.foodres.2013.10.007>
- Zhou L, Wang L, Palais R, Pryor R & Wittwer C T (2005). High-resolution DNA melting analysis for simultaneous mutation scanning and genotyping in solution. *Clinical Chemistry* 51(10): 1770-1777 <https://doi.org/10.1373/clinchem.2005.054924>



© 2022 by the author(s). Published by Ankara University, Faculty of Agriculture, Ankara, Turkey. This is an Open Access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.